

# CLOUD NINE Face and Body Sculpting Device

CLOUD NINE are fiercely kind, we celebrate beauty in all its forms; it is as individual as you are.

We provide the styling tools for you to express your most authentic self, they are designed with kindness in mind. Kinder to hair, to the environment, and to you.

For stronger, healthier hair that's uniquely you.

CLOUD NINE – BE FIERCELY YOU

## **How to use:**

Remove the protective cover from the device and switch on

There are four different modes for the device, and they are in the following order - Heating and Red light, Heating and Vibration, Cooling and Blue Light and Vibration and Blue Light. Press the button to manually change through each setting.

Each setting will select for 5 minutes. After 5 minutes the device will pause for 1 minute with a red light flash. Selecting the mode again will continue use, otherwise, the device will turn off.

Gently glide the device in an upward, outward motion to drain towards the lymph nodes by the ears correctly

Move to the jawline. Glide the device upwards along the jawline towards the ears for drainage. Or, continue to work upwards along the jawline to massage any tension

Move to the neck area and direct the device up towards the jaw starting from the collarbone

The device can also be used on the arms, legs, and other areas of the body to promote drainage and sculpting – always glide using a upward motion to aid drainage

To switch off press and hold the on/off button. The light on the device will go out when the device is off

## **Settings**

*Heating - helps to open the pores and increase the penetration and absorption of skincare products such as serums and moisturisers*

*Cooling – soothes and calms the skin. Tightens pores and reduces puffiness and inflammation.*

*Infra-red light - Infra Red Light emits beneficial light rays that energize cells and stimulate the skin's natural renewal process. Naturally plumping the skin by increasing collagen production and working to erasing fine lines and wrinkles*

*Ultrasonic vibrations – stimulates muscles and improves facial contouring*

*Ultrasonic Vibrations work to reduce inflammation and puffiness. Increasing muscle stimulation to improve facial contour by defining cheekbones and jawlines. Whilst also helping to reduce fine lines and improve skin texture.*

*Blue light – fights acne and blemishes*

*Blue light Therapy works to heal existing blemishes and prevent future breakouts. Improving skin texture, by destroying acne-causing bacteria, reducing pores and locking in moisture.*